

## Test Case Planning- FOR UC03- SPRINT 1

<b>Group members</b>	<b>Melinda Vo, Aidan, Ben, Benhi, Hunny, Kevin</b>
----------------------	--

### 1. Use Case 1 - <<name of the use case>>

<b>Test Case ID</b>	<b>TC01 - Correct comparison between healthy and user's data</b>
<b>Test Objective</b>	<b>Verify Healthy statistics are correctly compared</b>
<b>Preconditions</b>	<b>1. Login to the account 2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional)) 3. Connect to the watch</b>
<b>Test steps</b>	<b>1. Open the app 2. Click on Login button (if you don't have an account then press create new account) 3. Connect your watch to the application/app 4. Go to personal data 5. Tap on statistic</b>
<b>Input Values</b>	<b>1. Name: Lacy Lui 2. SSN: 123-123-1234 3. Physician: Doctor Candy</b>  <b>Data should be collected with watch</b>
<b>Expected results</b>	<b>If in healthy range:  **GRAPH OF USER'S DATA**</b>

	<p><b>Your heart rate/blood oxygen/ hours of sleep/ blood pressure is good</b></p> <p><b>If in unhealthy range:</b></p> <p><b>**GRAPH OF USER'S DATA**</b></p> <p><b>Your heart rate/blood oxygen/ hours of sleep/ blood pressure is bad</b></p>
--	--

<b>Test Case ID</b>	<b>TC02 - Display Graph</b>
<b>Test Objective</b>	<b>Verify User's data is correctly graphed</b>
<b>Preconditions</b>	<ol style="list-style-type: none"> <li><b>1. Login to the account</b></li> <li><b>2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional))</b></li> <li><b>3. Connect to the watch</b></li> </ol>
<b>Test steps</b>	<ol style="list-style-type: none"> <li><b>1. Open the app</b></li> <li><b>2. Click on Login button (if you don't have an account then press create new account)</b></li> <li><b>3. Connect your watch to the application/app</b></li> <li><b>4. Go to personal data</b></li> <li><b>5. Tap on statistic</b></li> </ol>
<b>Input Values</b>	<p><b>Name: Hunny Biguvu</b></p> <p><b>SSN:- 021-***-***9 (Optional)</b></p>

	<b>Physician: Doctor Candy(optional)</b>
<b>Expected results</b>	<b>Graph on blood pressure, blood oxygen, Heart rate, and hours of sleep</b>

<b>Test Case ID</b>	<b>TC03 - Verify values</b>
<b>Test Objective</b>	<b>Verify that the user's data does not equal to 0</b>
<b>Preconditions</b>	<ol style="list-style-type: none"> <li><b>1. Login to the account</b></li> <li><b>2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional))</b></li> <li><b>3. Connect to the watch</b></li> </ol>
<b>Test steps</b>	<ol style="list-style-type: none"> <li><b>1. Open the app</b></li> <li><b>2. Login with valid username + password</b></li> <li><b>3. Open Data tab</b></li> </ol>
<b>Input Values</b>	<b>Name: Lacy Lui</b> <b>SSN: 123-123-1234</b> <b>Physician: Doctor Candy</b>
<b>Expected results</b>	<b>Data set of blood pressure, heart rate, and blood oxygen.</b>

<b>Test Case ID</b>	<b>TC04 - User notification for Health Alerts</b>
<b>Test Objective</b>	- Verify the app is notifying appropriate notifications or alerts based on health history
<b>Preconditions</b>	<ol style="list-style-type: none"> <li>1. Login into account</li> <li>2. Enter your data</li> <li>3. Enable notifications</li> </ol>
<b>Test steps</b>	<ol style="list-style-type: none"> <li>1. Open the app</li> <li>2. Login with a valid username</li> <li>3. Open your profile</li> <li>4. Turn on the notifications for alerts</li> </ol>
<b>Input Values</b>	<ol style="list-style-type: none"> <li>4. Name: Lacy Lui</li> <li>5. SSN: 123-123-1234</li> <li>6. Physician: Doctor Candy</li> </ol> <p>Data should be collected with watch</p>
<b>Expected results</b>	<ul style="list-style-type: none"> <li>- Notifications should be sent as:</li> </ul> <p><b>!! Abnormal Heart Rate/Blood Pressure/Blood Oxygen/Hours of sleep !!</b></p>