

Test Case Planning

Group members	Melinda Vo, Aidan, Ben, Benhi, Hunny, Kevin
---------------	---

1. Use Case 1 - <<name of the use case>>

Test Case ID	TC01 - User Registration
Test Objective	Verify that a user can successfully register
Preconditions	<ol style="list-style-type: none">1. Valid email address and password2. Password must have at least 8 characters, a combination of uppercase, lowercase, numbers, and a symbol
Test steps	<ol style="list-style-type: none">1. Click on the registration page.3. Enter a valid email address4. Enter a valid password5. Confirm password6. Click on the "Register" button.
Input Values	<ol style="list-style-type: none">1. Normal registration Email: johndoe@gmail.com Password: ExamplePass&123 Confirm Password: ExamplePass&1232. Invalid registration Email: johndoe@gmail.com Password: ExamplePass&123 Confirm Password: ExamplePass&123

Expected results	<ol style="list-style-type: none"> 1. Account successfully created <ol style="list-style-type: none"> a. Confirmation message showing successful registration b. User redirected to login page 2. Password isn't meeting the requirements <ol style="list-style-type: none"> a. Missing number or @ 3. Already have an Account <ol style="list-style-type: none"> a. Email address or SSN matched
-------------------------	---

Test Case ID	TC02 - Display Graph
Test Objective	Verify User's data is correctly graphed
Preconditions	<ol style="list-style-type: none"> 1. Login to the account 2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional)) 3. Connect to the watch
Test steps	<ol style="list-style-type: none"> 1. Open the app 2. Click on Login button (if you don't have an account then press create new account) 3. Connect your watch to the application/app 4. Go to personal data 5. Tap on statistic
Input Values	<p>Name: Hunny Biguvu</p> <p>SSN:- 021-***-***9 (Optional)</p> <p>Physician: Doctor Candy(optional)</p>

Expected results	Graph on blood pressure, blood oxygen, Heart rate, and hours of sleep
-------------------------	--

Test Case ID	TC03 - Verify values
Test Objective	Verify that the user's data does not equal to 0
Preconditions	<ol style="list-style-type: none"> 1. Login to the account 2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional)) 3. Connect to the watch
Test steps	<ol style="list-style-type: none"> 1. Open the app 2. Login with valid username + password 3. Open Data tab
Input Values	Name: Lacy Lui SSN: 123-123-1234 Physician: Doctor Candy
Expected results	Data set of blood pressure, heart rate, and blood oxygen.

Test Case ID	TC04 - Account Deletion
---------------------	--------------------------------

Test Objective	<ul style="list-style-type: none"> - Delete account
Preconditions	<p>The User must have an account</p> <p>Must be logged in.</p>
Test steps	<ol style="list-style-type: none"> 1. Open App 2. Navigate to Account Settings 3. Tap Delete Account 4. Tap Confirm Deletion
Input Values	<ol style="list-style-type: none"> 1. Confirmation Input: Yes
Expected results	<ul style="list-style-type: none"> - Account is deleted - The user receives confirmation message (pop up) - The user is redirected to account sign-in page (with their login now deleted)