

## Test Case Planning

Group members	Melinda Vo, Aidan, Ben, Benhi, Hunny, Kevin
---------------	---

### 1. Use Case 1 - *Patients Access, tracking, and editing of data*

Test Case ID	TC01 - User Registration
Test Objective	Verify that a user can successfully register
Preconditions	7. Valid email address and password 8. Password must have at least 8 characters, a combination of uppercase, lowercase, numbers, and a symbol
Test steps	1. Click on the registration page. 9. Enter a valid email address 10. Enter a valid password 11. Confirm password 12. Click on the "Register" button.
Input Values	3. Normal registration Email: <a href="mailto:johndoe@gmail.com">johndoe@gmail.com</a> Password: ExamplePass&123 Confirm Password: ExamplePass&123 4. Invalid registration Email: <a href="mailto:johndoe@gmail.com">johndoe@gmail.com</a> Password: ExamplePass&123 Confirm Password: ExamplePass&123

<b>Expected results</b>	<b>4. Account successfully created</b> <ol style="list-style-type: none"> <li>a. Confirmation message showing successful registration</li> <li>b. User redirected to login page</li> </ol> <b>5. Password isn't meeting the requirements</b> <ol style="list-style-type: none"> <li>a. Missing number or @</li> </ol> <b>6. Already have an Account</b> <ol style="list-style-type: none"> <li>a. Email address or SSN matched</li> </ol>
-------------------------	---

<b>Test Case ID</b>	<b>TC02 - Display Graph</b>
<b>Test Objective</b>	<b>Verify User's data is correctly graphed</b>
<b>Preconditions</b>	<ol style="list-style-type: none"> <li><b>1. Login to the account</b></li> <li><b>2. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional))</b></li> <li><b>3. Connect to the watch</b></li> </ol>
<b>Test steps</b>	<ol style="list-style-type: none"> <li><b>1. Open the app</b></li> <li><b>2. Click on Login button (if you don't have an account then press create new account)</b></li> <li><b>3. Connect your watch to the application/app</b></li> <li><b>4. Go to personal data</b></li> <li><b>5. Tap on statistic</b></li> </ol>
<b>Input Values</b>	<b>Name: Hunny Biguvu</b> <b>SSN:- 021-***-***9 (Optional)</b> <b>Physician: Doctor Candy(optional)</b>

<b>Expected results</b>	<b>Graph on blood pressure, blood oxygen, Heart rate, and hours of sleep</b>
-------------------------	--

<b>Test Case ID</b>	<b>TC03 - Verify values</b>
<b>Test Objective</b>	<b>Verify that the user's data does not equal to 0</b>
<b>Preconditions</b>	<b>2. Login to the account</b> <b>3. Enter your data(Legal Name, SSN(optional), PHYSICIAN(optional))</b> <b>4. Connect to the watch</b>
<b>Test steps</b>	<b>1. Open the app</b> <b>2. Login with valid username + password</b> <b>3. Open Data tab</b>
<b>Input Values</b>	<b>Name: Lacy Lui</b> <b>SSN: 123-123-1234</b> <b>Physician: Doctor Candy</b>
<b>Expected results</b>	<b>Data set of blood pressure, heart rate, and blood oxygen.</b>

<b>Test Case ID</b>	<b>TC04 - Account Deletion</b>
---------------------	--------------------------------

<b>Test Objective</b>	<ul style="list-style-type: none"> <li>- Delete account</li> </ul>
<b>Preconditions</b>	<p>The User must have an account</p> <p>Must be logged in.</p>
<b>Test steps</b>	<ol style="list-style-type: none"> <li>1. Open App</li> <li>2. Navigate to Account Settings</li> <li>3. Tap Delete Account</li> <li>5. Tap Confirm Deletion</li> </ol>
<b>Input Values</b>	<ol style="list-style-type: none"> <li>2. Confirmation Input: Yes</li> </ol>
<b>Expected results</b>	<ul style="list-style-type: none"> <li>- Account is deleted</li> <li>- The user receives confirmation message (pop up)</li> <li>- The user is redirected to account sign-in page (with their login now deleted)</li> </ul>